

Live Viewer: Installation

Installation

- Download the live-viewer from the git repository: [Live Viewer in stash](#) and unzip it
- install python-zmq
- install python-inotifyx (or on fedora: pip install inotifyx (this needs python-devel to be installed))

(tested versions:

pyzmq version 14.0.1

libzmq version 4.0.4

inotifyx version 0.2.0)

Configuration

The live viewer as a configuration file for the sender and for the receiver located in **conf/**

Sender

Modify **sender.conf**:

- specify the monitored folder and the folder the "mv" should copy the data into:

sender.conf

```
# folder you want to monitor for changes
# inside this folder only the subdirectories "commissioning", "current" and
"local" are monitored
watchFolder      = ...
# Target to move the files into
cleanerTargetPath = ...
```

- specify the file format to be monitored:

sender.conf

```
# the formats to be monitored, files in an other format will be be
neglected
monitoredFormats = [".tif", ".cbf"]
```

- specify the hosts which are allowed to connect to the sender: e.g.

sender.conf

```
# list of hosts allowed to connect to the sender
receiverWhiteList = ["zitpcx19282", "zitpcx22614", "lsdma-lab04"]
```

Receiver for LiveViewer

Modify **receiverLiveViewer.conf**:

- specify the target folder for the ring buffer:

receiverLiveViewer.conf

```
# where incoming data will be stored to"
targetDir          = ...
```

- specify the ip where the sender is running: e.g.

receiverLiveViewer.conf

```
# IP of dataStream-socket to send signals back to the sender
senderComIp        = lsdma-lab04.desy.de
```

- specify the size of the ring buffer: e.g.

receiverLiveViewer.conf

```
# size of the ring buffer for the live viewer
maxRingBufferSize = 10
```